

Laser Handpiece and Core Module Portfolio

Integrated OEM Laser Delivery Solutions for Aesthetic Systems



Precision. Proven



About Us

Lasetica (Jinan) Co., Ltd. (established 2025) is a specialist manufacturer founded by a seasoned team from the handheld aesthetic laser industry, focusing on high-quality,rapid-response OEM/ODM services. Leveraging years of industry experience and deep insights into foundational technical principles, our core team provides global clients with highly customized solutions that far exceed industry standards, ensuring both exceptional quality and rapid delivery.

Technical Capabilities

Deep Expertise in Laser Engineering:

- 20+ years in Nd:YAG for aesthetic and dermatologic applications
- Expertise in active/passive Q-switching and frequency doubling
- Energy stability ≤5% RMS

Product Portfolio

Versatile Handpiece Solutions:

- N1 Standard passive Q-switch; costeffective
- D1 High-Energy passive Q-switch; for deep pigment removal and large-area treatments
- S1 Precision EO (active) Q-switch; precise control

Section 2 Customization Capabilities

OEM/ODM partner enabling rapid innovation.

- Full customization: parameters and appearance
- Accelerated prototyping: technical proposal and quote within 10 days; samples in 8–12 weeks
- Your IP remains confidential under NDA

Global Support

Global reach, dependable support.

- 30+ countries served: a network you can rely on
- Rapid response: spare parts dispatched within 48 hours; technical support 10 hours/day, 7 days/week
- Team empowerment: training and documentation included



Leo
Laser R&D Expert (20+ Years)
Leads core technology
architecture for precision and
energy stability.



Ryan
Operations & Supply Chain
Lead (18+ Years) ·
Engineering a customercentric network for rapid
delivery and reliable service.



Duncan
Senior Industrial Engineer
Fusing ergonomic design
with process optimization for
long-term reliability.

Customization Services

Every customer's needs are unique. Lasetica delivers both standardized products and **customized solutions**, from parameter adjustment to full OEM/co-development, tailored to your specific requirements.

Range Parameter Customization

Customize core parameters—energy output, pulse width range, energy density, specific wavelengths—to precisely match end-user treatment needs.

Optical Customization

Customizable spot size, spot uniformity (RSD), wavelength combination to meet diverse treatment scenarios.

♣ Deep OEM Cooperation

Complete solution design, exclusive technology licensing, end-to-end customization.

Industrial Design & Appearance

Logo engraving, color customization, integrated internal structure and enclosure design—create exclusive products that match your brand image.

Description Interfaces

Flexible adaptation to integrated systems electrical interfaces, control protocols, cooling interfaces and other technical specifications.

♦ Fast Delivery

Standard: 2–3 weeks; Customized: 4–6 weeks; Deep OEM: 8–12 weeks; NDA supported

Specification Comparison

Parameter	N1 (Standard)	D1 (High Energy)	S1 (Precision)
Wavelength	1064/532 nm	1064/532 nm	1064/532 nm
Energy Output	180-230 mJ (single pulse)	≤450 mJ (single pulse)	≤450 mJ
Pulse Width	8–12 ns	4.5–12 ns	6–10 ns
Peak Power (estimated)	12–25 MW	20–40 MW	25–50 MW
Energy Stability	≤5% RMS	≤5% RMS	≤4% RMS
Cooling Flow	2-3.6 L/min	2–3.6 L/min	2-3.6 L/min
Recommended Applications	Entry-level equipment, standard pigment treatment	Deep pigment removal, large- area treatment	advanced treatments, research applications

Optional Accessories

Accessory Name	Specification	Compatible Models
532 nm Frequency Doubling Head	Green light output, pigment treatment enhancement	N1/D1/S1
Replaceable Focus Head Kit	3 spot sizes (2/3/5 mm)	N1/D1/S1
Honeycomb Spot Treatment Head	RSD ≤10%, uniformity enhancement	N1/D1/S1
Pixel Spot Treatment Head	7×7 pixel spot array	N1/D1/S1
External Trigger Interface	TTL signal, research applications	S1 only

S1 Series

Precision Control Flagship









Key Product Features

Ultra-Narrow Pulse Width

6–10 ns ultra-narrow pulse, EO (active) Q-switch technology delivers selective photothermolysis for targeted interaction while minimizing collateral thermal impact.

Full Parameter Adjustability

Energy and frequency independently adjustable—energy 10–450 mJ, frequency 1–15 Hz, pulse width 6–10 ns—adapting to various treatment protocols.

Superior Stability

≤4% RMS energy stability; 8-hour continuous operation fluctuation ≤3%, supporting consistent output and system reliability.

Premium Applications

Ideal for high-end aesthetic equipment, research applications, precision treatment scenarios—meeting strict requirements for energy and frequency control.

Product Parameters

 Wavelength
 1064 nm (532 nm)
 Energy Output
 ≤450 mJ

 Pulse Width
 6–10 ns
 Peak Power
 25–50 MW

 Energy Stability
 ≤4% RMS
 Cooling Flow
 2–3.6 L/min

Note: Unless otherwise stated, parameters are measured at 25 ± 2 °C, rated cooling, standard optical interface conditions; maximum energy typically corresponds to low-frequency mode; peak power is an estimated value.

D1 Series

High-Energy Platform for Deep Treatments







Key Product Features

High Energy Output

Single pulse energy up to 450 mJ, providing sufficient energy margin for deep pigment removal and large-area treatment.

Proven Reliability

≤5% RMS energy stability, validated by factory burn-in testing (48–72 hours) for long-term stable operation.

Large Spot Coverage

Equipped with 6–8 mm large spot focus head, adapting to various treatment scenarios, significantly improving work efficiency.

Professional Applications

Ideal for deep pigment removal, large-area treatments, professional aesthetic systems requiring high energy.



Product Parameters

 Wavelength
 1064 nm (532 nm)
 Energy Output
 ≤450 mJ

 Pulse Width
 4.5–12 ns
 Peak Power
 20–40 MW

 Energy Stability
 ≤5% RMS
 Cooling Flow
 2–3.6 L/min

Note: Unless otherwise stated, parameters are measured at 25 ± 2 °C, rated cooling, standard optical interface conditions; maximum energy typically corresponds to low-frequency mode; peak power is an estimated value.

N1 Series

Cost-Effective Standard Platform







Key Product Features

Efficient Energy Output

Peak single-pulse energy up to 230 mJ, covering most standard pigment treatments.

Fast Delivery

Standardized platform and mature processes shorten the lead time from technical proposal to mass production.

Easy Maintenance

Modular design for convenient maintenance, reducing equipment downtime.

Recommended Applications

Entry-level aesthetic systems; standard pigment treatments; costeffective market positioning.



Product Parameters

 Wavelength
 1064 nm (532 nm)
 Energy Output
 180–230 mJ

 Pulse Width
 10–12 ns
 Peak Power
 12–25 MW

 Energy Stability
 ≤5% RMS
 Cooling Flow
 2–3.6 L/min

Note: Unless otherwise stated, parameters are measured at 25 ± 2 °C, rated cooling, standard optical interface conditions; maximum energy typically corresponds to low-frequency mode; peak power is an estimated value.

Customization Service Plans

Standard Version

Delivery time: 2-3 weeks

- · Direct supply of stock models
- · Fast delivery
- · Standard warranty

Customized Version

Delivery time: 4-6 weeks

- Parameter fine-tuning
- Appearance customization
- Logo engraving
- · Packaging customization

Deep OEM Customization

Delivery time: 8-12 weeks

- · Complete solution design
- Exclusive technology licensing
- · End-to-end customization
- · Production line support

Cooperation Process

1

Requirement Communicati on

Sign NDA

② Solution Design

10-day technical solution + quotation

③ Sample Verification

2–3 weeks sample delivery On-site testing support **4 Mass Production**

Small-batch production from 3 weeks Rolling delivery

After-Sales Support System

Warranty

1 year or 1 million shots Whichever comes first.

Response

Emergency failures Response within 4 hours

Spare Parts

Mainland China Dispatched within 24–48 hours

Support

Technical support 7 days/week Remote diagnosis

Training

Operation training videos Maintenance guides

Upgrade

Free firmware updates Upgrade services

Contact Us

Address: Room 1-711, Xinyuan Guoji, No. 36 Gongyenanlu, Jinan, Shandong 250000, China Website: www.lasetica.com | Business: info@lasetica.com | Technical: tech@lasetica.com



Website



Email



Legal Statement and Disclaimer

Disclaimer: Technical parameters and product descriptions in this catalog are for reference only; actual specifications are subject to the final contract. Lasetica reserves the right to modify product specifications without prior notice.

Usage Restrictions: This catalog content is protected by copyright and may not be reproduced, distributed, or used for commercial purposes without written permission from Lasetica.

Liability Limitation: Lasetica is not responsible for direct or indirect losses caused by use of information in this catalog. Product applications must comply with local laws, regulations, and medical device requirements.

Intellectual Property: Trademarks, logos, and product names in this catalog are the intellectual property of Lasetica or its affiliated companies, protected by Chinese and international intellectual property laws.

Certification Statement: Product certification status is subject to official certification body certificates. Expressions such as "CE in preparation" and "ISO 13485 in preparation" represent current progress only and do not constitute certification commitment.

Governing Jurisdiction: This catalog is governed by the laws of the People's Republic of China; any disputes shall be under the jurisdiction of the People's Court of Jinan City. For legal or compliance issues, contact: legal@lasetica.com

嗲 Lasetica | 来光科技

Your Partner in OEM/ODM Laser Handpiece Solutions

Your Partner for Life

Address: Room 1-711, Xinyuan Guoji, No. 36 Gongyenanlu, Jinan, Shandong 250000, China

Website: www.lasetica.com | Business: info@lasetica.com | Technical: tech@lasetica.com





Wehsite

Email

Copyright © 2025 Lasetica (Jinan) Co., Ltd. | Last Updated: 2025-10 | www.lasetica.com